

# **900 SERIES DELIVERY CART SERVICE MANUAL**

**INCLUDES: OPERATION  
MAINTENANCE  
TROUBLESHOOTING  
INSTALLATION**

**VERSION 1  
DATED: AUGUST 1, 2003**



**SUMMIT DENTAL SYSTEMS**

**5340 NW 10<sup>TH</sup> TERRACE  
FT. LAUDERDALE, FL 33309  
(800) 275-3368  
(954) 938-6831 \* FAX (954) 938-6833  
[WWW.SUMMITDENTAL.COM](http://WWW.SUMMITDENTAL.COM)**

# Table of Contents

<b>WARRANTY.....</b>	<b>3</b>
<b>IMPORTANT INFORMATION.....</b>	<b>4</b>
<b>Weights &amp; Dimensions.....</b>	<b>5</b>
<b>Operator's Instructions.....</b>	<b>6</b>
<b>Flow Diagram.....</b>	<b>7</b>
<b>PURGING YOUR SYSTEM.....</b>	<b>8</b>
<b>TROUBLESHOOTING GUIDE.....</b>	<b>9</b>
<b>INSTALLATION INSTRUCTIONS.....</b>	<b>10</b>
<b>Utility Center Plumbing Layout.....</b>	<b>11</b>
<b>Utility Center Tubing Diagram.....</b>	<b>12</b>
<b>Cabinet Parts List.....</b>	<b>13</b>

## WARRANTY

Summit Dental Systems (SDS) warrants its products against defects in materials or workmanship from the date of shipment to the Buyer as follows:

**Summit Dental Systems (SDS) Equipment:**

Handpiece Control Assembly, Utility Center  
Cabinet, Removable Plastic Counter Tray

**Warranty Period:**

3 Years  
Date of Shipment  
F.O.B. Ft. Lauderdale, Florida

Summit Dental Systems' sole obligation under the warranty is to provide parts for repair, or at its option, to provide a replacement product (excluding all labor and shipping fees). The buyer shall have no other remedy. All special, incidental, coincidental damages and normal wear and tear to the cabinets are excluded. The warranty does not cover damage from improper installation or maintenance, accident or misuse. The warranty does not cover damage resulting from the use of cleaning disinfecting or sterilization chemicals and processes. Failure to follow instructions provided in Summit Dental Systems' Operation and Installation Manuals (Owner's Guides) may void the warranty. In the event Warranty service must be performed to correct any defect, only an authorized Summit Dental System dealer may perform any and all warranty repairs. Any repairs by unauthorized dealers, technicians, or repairmen may void the warranty.

- In the case of a defective warranty item, a copy of the replacement invoice, serial number of the product under which it was replaced, and a description of symptoms of the defect must be returned with the part to Summit Dental Systems, 999 NW 53<sup>rd</sup> Court, Ft. Lauderdale, FL 33309, USA, in order to receive credit. Any and all expenses for freight, labor to perform warranty service, and purchase of spare parts is the responsibility of the buyer. Any fraudulent claims made may void the warranty. Any additional warranty that may be provided by an authorized Summit Dental Systems dealer is the sole responsibility of said dealer.

- A 15% restocking fee will be applied to the authorized dealer's account for any parts returned do to a misdiagnosis by the service technician performing the warranty service.

- SDS reserves the right to make changes or improvements on any products without being required to modify existing equipment in a like manner.

**Please complete and retain for your records the following Information:**

In case of warranty part replacement/repair or when ordering a part, please call your authorized Summit Dental Systems dealer and have the following information available:

<b>Owners' Name:</b>	<b>Phone #:</b>
<b>Model #:</b>	<b>SDS Serial #:</b>
<b>Dealer:</b>	<b>Phone:</b>
	<b>Purchase Date:</b>

# Important Information

---

## Technical Specifications

**Air Regulator Pressure: 80 PSI**

**Water Regulator Pressure: 40 PSI**

**Handpiece Pressure: 30-35 PSI (for most high-speed handpieces)  
40-45 PSI (for most low-speed handpieces)**

**CAUTION:** When adjusting the handpiece pressure, do not over tighten the screws. This may result in unnecessary damage to the handpiece control block.

- Pay careful attention when unpacking the delivery cart and its accessories. *Damage caused by mishandling the equipment during unpacking or installation is not covered under warranty.*
- **New Owner:** Please read, sign and submit the warranty registration form that is located on page 3 of this manual. Failure to return this form may void the warranty.
- Read this manual carefully. If you have any questions, please call **SDS Technical Service** at **1-800-275-3368**.

### Product Label Location


**The product label can be located inside the handpiece control door.**

<b>SUMMIT DENTAL SYSTEMS</b> Ft, Lauderdale, FL USA		
MODEL No: _____		
SERIAL No: _____		
PRODUCT: _____		
VOLTS		HERTZ
		AMPS

# Weight & Dimensions

---


## Shipping Information



*Model 900 Series Delivery Carts*

Shipping Weight: 130 lb.


## General Dimensions



Cabinet Weight: 88 lb.

# Operator's Instructions

---



## A. MASTER ON/OFF SWITCH

Opens and closes the valves in the floor-mounted Utility Box, supplying air to the complete system. Shut off nightly.

## B. WATER ON/OFF SWITCH

Opens and closes the valves in the floor-mounted Utility Box, supplying water to the complete system. Shut off nightly.


## C. WATER COOLANT FLOW CONTROL

Controls the flow of water coolant to handpieces. Turn counter-clockwise to increase, and clockwise to decrease or to fully cease for low-speeds.

## D. HANDPIECE PRESSURE GAUGE

Indicates pressure adjusted by the Drive-Air Control Knob. Avoid exceeding the manufacturer's recommended handpiece pressure. (30 – 35 PSI for most Highspeed Handpieces)  
(40 – 45 PSI for most Low-speed Handpieces).

# 900 Series Flow Diagram



# Purging Your System

---

## Daily Maintenance

**Purging with air:** At the end of every day, the lines should be purged with air to prevent the growth of biofilm.

1. Remove handpieces from tubing.
2. Empty the water bottle, then reinstall it (if water bottle is applicable, when using city water, turn off water from the source).
3. Hold the handpiece tubing and syringe over a pail. Turn the unit on, wait a few moments and then operate the flush toggle (if applicable), syringe and foot control until water is purged from the system.
4. Turn unit off.

**Disinfecting the bottle (if applicable):** Fill the bottle with 100 ml disinfectant solution, shake vigorously and let settle for 10 minutes. Shake again, and then rinse

## Weekly Maintenance

The weekly cleaning procedure should be performed at least once a week, preferably at the start of the week before treating patients. If the unit is stored for any length of time, perform the weekly maintenance routine immediately before and after storage.

1. Purge the unit with air (see Daily Maintenance).
2. Flush the system with disinfectant solution:
  - a. Turn unit off. Empty the water bottle (if water bottle is applicable), replacing the water with cleaning solution (see Disinfectant Solution).
  - b. Remove handpieces from tubing and hold the handpiece tubing and syringe over a pail. Turn unit on, wait a few moments and then operate the flush toggle (if applicable), syringe and foot control until a continuous stream of solution is running through the system.
3. Allow the disinfectant to remain in the unit for at least 10 to 20 minutes and then flush the system again until all the solution is used up.
4. Purge the unit with air:
  - a. Hold the handpiece tubing and syringe over a pail. Turn the unit on, wait a few moments and then operate the flush toggle (if applicable), syringe and foot control until all solution is purged from the system.
  - b. Turn unit off. (If the unit is to be stored, stop here.)
5. Fill with clean water:
  - a. With the unit turned off, remove the empty disinfectant bottle (if applicable). Replace with clean bottle and water.
  - b. Hold Remove handpieces from tubing and hold the handpiece tubing and syringe over a pail. Turn unit on, wait a few moments and then operate the flush toggle (if applicable), syringe and foot control until a continuous stream of solution is running through the system.

**Disinfectant Solution:** Use 100 ml (9 parts tap water & 1 part 5.25% Sodium hypochlorite (household bleach) of disinfectant solution for each application per week. Always use a fresh mixture every week.

# Troubleshooting Guide

---

## SYMPTOM

1. Handpiece Lacks Power

## PROBABLE CAUSE

- A. Check regulator adjustment (80 psi).
- B. Check handpiece pressure adjustment on control block.
- C. Plugged air filter.
- D. Pinched supply tubing. Check for kinks.
- E. Bad handpiece gasket at connection with tubing.
- F. Defective handpiece.

2. Water coolant does not shut off when Foot Control is released.

- A. Adjust air pressure to 80 psi – water pressure to 40 psi.
- B. Foot Control is not exhausting.
- C. Defective water relay in valve.

3. More than one handpiece is operating.

- A. Handpiece is not completely in the hanger.
- B. Improper adjustment of Pilot Valve in the hanger.
- C. Kinked or pinched signal line from the Pilot Valve.

4. Insufficient water coolant.

- A. Adjust coolant flow valve.
- B. Water filter may be plugged.
- C. Plugged handpiece.
- D. Kinked or pinched tubing.
- E. Improper adjustment of water relay.

5. Water coolant is running from handpiece while in holder.

- A. Water pressure is too high.
- B. Air pressure is too low.
- C. Handpiece holder out of adjustment.



6. Water coolant is running continuously.

- A. Purge switch is on (on applicable systems).
- B. Water pressure is too high.
- C. Handpiece holder out of adjustment.
- D. Improper adjustment of water relay.

*If you continue to experience difficulties, please call Technical Service Team at 1-800-275-3368.*

# Installation Instructions


---



1. Your 900 Series Delivery Carts are pre-assembled for your convenience. Just open the box, take out the cart, place it on the floor and hook-up the utility center (see pages 11& 12).

# Utility Center Plumbing Layout


(Not to Scale)




NOTE: All pipes must be terminated 1 inch above floor.  
 Position box so that 2" opening faces post side.

# Utility Center Tubing Diagram

---




*figure 1*  
*Air side*



*figure 2*  
*Water Side*

## Cabinet Parts List

---



Description	Part Number
A. Removable Plastic Surface Tray	1-020-0101
B. Aluminum Handles	5-020-0034
C. Breadboard	2-020-XXXX ( <i>specify color</i> )
D. Small Storage Drawer	2-020-XXXX ( <i>specify color</i> )
E. Medium Storage Drawer	2-020-XXXX ( <i>specify color</i> )
F. Large Storage Drawer	2-020-XXXX ( <i>specify color</i> )
G. Dual Caster Wheels	2-020-0106
H. Cart Frame	2-020-0020
I. Handpiece Control Door	2-020-XXXX ( <i>specify color</i> )
J. Door Hinge ( <i>not shown</i> )	5-020-0020
K. Drawer Glide Kit ( <i>not shown</i> )	2-020-0108