



*Figure 25*

**Backrest Parts  
Model - 3050**

## Backrest Parts List - Model 3050

*refer to Figure 25 - page 43*

<b>Item</b>	<b>Part #</b>	<b>Description</b>
10	4-070-0029	Washer, Smooth, 5/32"
25	2-010-0055	Nut, 5/32"
43	2-010-0056	Pin, Operating Support
60	2-010-0170	Support, Rubber, Backrest, 3050/1750
61	2-010-0171	Spring, Pressure Plate, 3050/1750
67	4-010-0010	Bolt, Pressure Plate, 3050/1750
83	2-010-0172	Plate, Headrest, Pressure 3050/1750
84	2-010-0173	Frame, Backrest, 3050/1750
85	2-010-0174	Support, Headrest, Pressure, 3050/1750
86	2-010-0175	Frame, Headrest, Pressure, 3050/1750
87	4-010-0011	Bolt, Articulation, Backrest
88	2-010-0176	Headrest Guide Assembly, 3050/1750