



Figure 32

**Foot Control Parts
Model 1750**

Foot Control Parts List - Model 1750

refer to Figure 32 - page 59

Item	Part #	Description
01	2-010-0076	Frame, Foot Control
02	2-010-0077	Cover, Foot Control
03	2-010-0057	Lever, Foot Control
04	2-010-0079	Collar, Control Lever, FC,
05	7-010-0004	Mount, Micro-Switches
06	2-010-0080	Collar, Control Lever, Retainer, FC
07	4-070-0038	Retainer, (C), FC,
08	7-010-0005	Button, Entry/Exit, FC,
09	2-010-0081	Spring, FC, Entry/Exit Buttons
10	2-010-0082	Retainer, (C), FC, Entry/Exit Button
12	4-070-0040	Screw, Phillips, SM, 1/8" X 1/2"
13	7-010-0002	Switch, Micro