



Figure 26

**Headrest Parts
Model - 3050**

Headrest Parts List - Model 3050

refer to Figure 26 - page 45

| Item | Part # | Description |
|-------------|---------------|--|
| 01 | 2-010-xxxx | Cushion, Headrest |
| 02 | 2-010-0340 | Double Articulating Headrest Assembly |
| 03 | 4-010-0049 | Screw, Phillips Pan Head 10-24 X 3/8" |
| 04 | 2-010-0341 | Cover, Headrest, Double Articulating |
| 05 | 4-010-0050 | Screw, Phillips Pan Head 10-24 X 5/16" |