

# BISCAYNE E.L.

## CHAIR



The Biscayne E.L. Chair delivers quality and value ideal for hygienists and dentists. With a clean look and simple design, the Biscayne E.L. is practical for any dental setting. Economically priced without sacrificing quality, the Biscayne E.L. is a viable option for dentists who want a functional chair they can rely on for their practice.

As with all SDS chairs, the electro-mechanic chair is covered under our industry leading 5-year warranty.

## FEATURES

- Steel base & backrest assembly
- Powder-coated baked finish
- Electro-mechanic base & back
- 1 programmable pre-position
- Automatic return
- Foot control
- Double articulating headrest
- Armrests
- 1 Outward rotating armrest

## SPECIFICATIONS

RISING CAPACITY	440/200kg
NET WEIGHT	230 lb (104 kg)
CHAIR WIDTH	26.5" (673mm)
CHAIR LENGTH W/ BACKREST DOWN	77" (1956mm)
MAX HEIGHT W/ BACKREST DOWN	34.5" (876mm)
MIN HEIGHT W/ BACKREST DOWN	20" (508mm)
MAX HEIGHT W/ BACKREST UP	57" (1448mm)
MIN HEIGHT W/ BACKREST UP	42" (1067mm)

# BISCAYNE E.L.

## HYGIENE LINE



The Biscayne E.L. Hygiene Package is the perfect fit for all your dental hygiene needs. Streamlined for success, the Hygiene Package comes equipped with the 1335 Biscayne E.L. Light, the 1402 2HP Automatic Delivery Unit, and either the 1615E Assistant's Arm or the 1607 Cuspidor (both with SE and HVE). With these standard features and additional optional features, you'll see why the Biscayne E.L. Hygiene Package is the economical solution to bring quality and functionality to your practice.

- A** Biscayne E.L. Chair
- B** 1335 Light w/ Post
- C** 1402 2HP Automatic Delivery System
- D** 1615 Economy Assistant's Arm w/ SE and HVE



**SDS**<sup>®</sup>  
SUMMIT DENTAL SYSTEMS

800.275.3368 | Tel: 954.730.3636 | Fax: 954.730.3602  
6277 NW 28th Way | Fort Lauderdale, FL 33309  
[www.summitdental.com](http://www.summitdental.com)